



CPAP BOX		PULMONDYNE O2 MAX	
KEEP SAME NUMBERED ITEMS IN THEIR RESPECTIVE BAGS  DON'T SHARE ITEMS WITH OTHERS. CONTACT INSTRUCTOR IF THERE ARE QUESTIONS	1	CPAP MASK WITH STRAPS (A)	ASSURE FOREHEAD SUPPORT IS CONNECTED
	1	ADAPTER (MAY BE BLUE) 22mmF X 30mmF (B)	IS UTILIZED WHEN YOU ARE USING NEBULIZER WITH CPAP. ATTACHES TO MASK
	1	"T" VALVE (C)	ATTACHES TO VALVE ON MASK
	1	ADJUSTABLE PEEP VALVE (D)	ATTACHES TO OPPOSITE END OF "T" VALVE. ADJUST ACCORDING TO EMS PROTOCOL
	1	CPAP WHITE CORRUGATED TUBE (E)	CAN BE LENGTHEN TO 72" ATTACHES TO PEEP VALVE
	1	NEBULIZER (F)	ATTACHES TO DOWN PORT OF "T" VALVE
	1	O2 SUPPLY TUBING (G)	CONNECT LITER FLOW FOR NEBULIZER
	1	FIXED FLOW GENERATOR (H)	ATTACHES TO DISS PORT
REVIEW EMS1 VIDEO ON YOUTUBE- <a href="#">INNOVATION ZONE-PULMONDYNE O2 MAX SYSTEM, AND PULMONDYNE SET UP</a>			